

Location Information in Numera EverThere

Most Recent Location

Picks up a GPS signal every 10 minutes*

Numera Libris mPERS device provides location information based on the satellite fix. The location information is uploaded into the EverThere cloud platform.

EverThere, consumer portal, is designed to provide device and user information to the caregivers. EverThere is not a location tracking tool.



A Satellite Fix is Precise

A satellite fix uses satellite information to detect the subscriber's location and is quite precise. It is denoted by a light blue transparent circle on the map and can provide a physical address. The Libris needs a clear view of the sky to pick up a satellite fix, so if the Libris is in a building, it is possible that a satellite fix cannot always be detected.

A Network Fix Can Cover a Wide Area

A network fix uses the closest cellphone tower to determine the subscriber's approximate location. It is denoted by a red circle on the map, can cover a wide area, and is less precise in providing the exact location.

Location Data Updates Every 10 Minutes

Under optimal conditions the Libris device picks up the location via satellite every 10 minutes. After every 12 events all the collected event data is automatically updated on EverThere, or it can be requested on demand by pressing "refresh".

*when Libris device has access to satellite or cellular network.

Location In a Case of an Emergency

In an emergency event, like in a case of a fall, the device will present either a satellite or network fix. In addition to the location information provided by the Libris device, the cellular-based 2-way cellular voice is the best way to determine the exact location of the patient.

How the Location Information is Displayed in EverThere

The EverThere service allows the user's care-team members to view the event activity list on the EverThere dashboard screen. The most current location of the Libris device is displayed on the map with a corresponding address. EverThere user can also see whether the location is a satellite or network fix. There is also an icon "refresh location" on the address screen which can be used to get the most recent location of the Libris device.